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
**DUST EXTRACTION FOR DENTAL LABORATORY
SYSTEMES D'ASPIRATION POUR PROTHESISTES
ABSAUGER FÜR DEN ZAHNTECHNIKER
ASPIRACIÓN PARA USO ODONTOTÉCNICO**

- MANUALE ISTRUZIONI
- OPERATOR'S HANDBOOK
- MANUEL D'UTILISATION
- GEBRAUCHSANWEISUNG
- MANUAL DE INSTRUCCIONES

CATTANI S.p.A.

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INDICE

Pagina

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DUST EXTRACTION FOR DENTAL LABORATORY

INDEX

Page

— INTRODUCTION	10
— WARNING AND SIGNALS	10
— STOCKS	10
— PROTECTIVE SHIELD	10
— ARM RESTS	11
— ASPI-LABOR AND MONO-LABOR	11
— MOTORIZED ELECTROVALVES	11
— FILTERS	11
— FILTER FUNCTIONING	11
— AVAILABLE FILTERS	12
— ASPIRATION ASSEMBLY	12
— AVAILABLE ASPIRATION UNITS	12
— CENTRALISED SYSTEM	12
— SAFETY IN WORK PLACES	13
— CONTROL PANELS AND TECHNICAL ROOM IN A DENTAL LABORATORY ...	13
— MAINTENANCE	14
— IMPORTANT NOTICE	14
— TRANSPORT AND STORAGE	15

Introduction

This is meant to illustrate the appliances, the systems, and their functioning to dental technicians and to installation engineers, to advise about maintenance and to inform engineers and users about hazards and necessary precautions for preventing accidents.

Signals

- **Electrical shock risk: even 230 V power can be lethal**

- **High temperature**

- **Generic danger sign**

- **Compulsory direction of flow and rotation sense**

Not all hazards and compulsory instructions can be indicated by signals, so the user must read all notices and consider them carefully. Failing to observe a signal or a notice can be harmful for the operator or damage the appliance. Protection devices must not be removed. Structural or functioning features of the appliance or of the system must not be modified.

Stocks (draw. page 40)

Two models of stock are available: one fits on a 30mm Ø pipe (fig. 1 det. A), the other on a 40mm Ø pipe (fig. 2 det. A). A 30 mm Ø pipe (fig. 1) is less cumbersome and it can be easily placed on any workbench. Through a 40mm Ø pipe (fig. 2) there is a larger passage for air and aspiration is more powerful, but often the central drawer of the bench must be removed. Both stocks are fitted with a filter to recover precious metal filings (fig. 1 and 2 det. F).

Protective shield (draw. page 41)

Protective shield can be very useful only if it is very flexible (fig. 3 det. A-B-C-D). It protects the operator from fragments that aspiration can not drawn in. Fig. 4 shows the shield fixed under the cover bench, where it slides off when not used.

Arm rests (draw. page 40-41)

Arm rests help the operator to get a correct working posture. Model A (fig. 5) can be used with the shell (fig. 2 det. B). Longer arm rest model B (fig. 5) facilitates working on stock (fig. 1 and 2 det. A)

Aspi-Labor and Mono-Labor (draw. page 42)

General Features: Aspi-Labor and Mono-Labor are two mobile aspirators for dental laboratory. They can easily moved from a working position to another; noiseless, they are manufactured for continuous service.

They are meant to aspirate resin and chalk dusts, metal filings etc. Precious metal filings can be easily recovered, allowing low costs. These appliances safeguard respiratory system of the operator.

Functioning: press on the main switch (connection to the power supply) and the unit one (fig. 6 det. 7) (I = ON/ O =OFF) to start aspiration. Gold filter must be put in place only to recover precious filings and must also be often emptied, to avoid a decrease in aspiration power.

When only dust is extracted, the dusty air passes through the stock or the arm and gets to the filter-bag where dust is trapped, while air is conveyed to the aspiration unit and then exhausted outside the appliance.

For more information, please refer to the specific operator manual about Aspi-Labor and Mono-Labor.

Motorised electrovalves (draw. page 43)

Motorised electrovalves make each working position independent from others connected to the same system. When the first operator pushes on the aspiration starting button, the electrovalve opens and switches on the first aspiration unit. The last operator stopping the aspiration will also switch off the last aspiration unit.

Electrical diagram of the system, page 43 fig. 7.

Filters (draw. page 43-44)

The different types of filter-bags (fig. 8) are not appreciated any more: dust extracted in a dental laboratory is mostly chalk dust, for which it not easy to find an efficient and at the same time long –lasting filter bag.

Chalk dust sticks to the walls of a still empty filter thus throttling the aspiration like a full filter does. Too frequent replacements of filter entail time and money spending which can weigh on laboratory budget.

To solve the problem, we have studied a pneumatic cleaning system of the filter (fig. 9 det. A) which is automatically carried out any time the aspiration stops. This way every time the aspiration is switched on, the filter is always clean. Pneumatic cleaning assures a endless life to filters.

Filter functioning (draw. page 44)

Through the aspiration pipe (fig. 9, det B) dusty air goes into the dust separator tangent to its wall, dust is then collected at the bottom of the container (E), whilst air rises from outlet (D), passes through the filter (A) and flows towards the aspiration

unit. Every time the aspiration is switched off, a timer opens the electrovalve (C) and a few second's air blow, with a pressure of about 6 bar, cleans the filter pneumatically.

Available filters (draw. page 45)

We manufacture different models of filters:

Item 090070 (fig. 10) to be coupled to aspiration units Uni-Jet 75, Tecno-Jet, and Flux-Jet; it can be used for stock aspiration outlet, for dry model polishing machine and retouch box.

Item 090080 (fig. 11), fitted with aspiration unit, especially recommended for polishing and sand blaster machines, and for retouch box.

Item 091050 and 091060 (fig. 12 and 13) to be used for middle and large-sized centralised systems, up to 24 working positions for each filter.

Aspiration assembly (draw. page 45)

The aspiration assembly is fitted with one or more aspirators plus a dust collecting filter (fig. 14).

Filter capacity, dimensions and number of aspirators are calculated according to the working positions, to stocks and to other possible aspiration outlets. The assembly is noiseless and can be installed in the working room.

Available aspiration units (draw. page 46)

Uni-Jet 75 for one working position - stock aspiration outlet - aspiration unit on the spot.

Tecno-Jet for one working position - stock aspiration outlet - aspiration unit at 8-10 m.

Flux-Jet for two stocks aspirating at the same time - aspiration unit at 8-10 m.

Medio-Jet for three stocks aspirating at the same time - aspiration unit at 10-12 m. max.

Maxi-Jet 2S for four stocks aspirating at the same time - aspiration unit at 10-15 m.

Uni-Jet 500 for six stocks aspirating at the same time - aspiration unit at 20 m.

Uni-Jet 1000 for twelve stocks aspirating at the same time - aspiration unit at 25 m.

Aspiration units can be completed with the available control panels, silencers and dampers.

Furthermore, smallest aspiration units can be fitted with box for sound-deadening and safety protection. Aspiration units can be used for continuous service, without stops.

It is forbidden closing the aspiration to head values above the ones stated for each aspiration unit and changing the amperage set on the thermal magnetic protector on the control panel for each aspiration unit or compressor head.

Centralised system (draw. page 45-47)

A centralised system allows the installation of the aspiration units and of the dust collecting tanks out of the working room, which is a special requirement of middle-or large-sized dental laboratories.

Fig. 15 shows a middle-sized laboratory with 10 aspiration outlets (C) on stocks aspiration piping with their diameters (A), the low tension control lines (D).

We recommend heavy PVC aspiration piping, glued following the directions of the pipe manufacturer. Piping can be fixed at the ceiling with descending pipes of 50

mm Ø every two or three working positions or to the walls, if there is room enough behind the bench drawers. As shown in the example, when placing the aspiration piping, it is important to follow carefully the instructions concerning the diameter progression, the gradient (45°) of all branch connections (B) to the main piping, and the tightness of piping: any throttling and leakage can cause noticeable frictions and flow reduction. To absorb harmful vibrations, benches and aspiration units must be connected to the system piping by means of flexible hose.

Motorised electrovalves can be installed on each working positions also when the filter is installed for way. When installing them, pay attention not to obstruct the holes on the fixing away of the electrovalve, respect the flow direction and place top of the valve up, as shown on the electrovalve itself. On these conditions, the aspirated dusty air, which overcomes the filter placed on the stock, can pass through the electrovalve.

The two aspiration units Uni-Jet 500 shown in this example are enough to service twelve positions working at the same time.

A centralised system is usually connected to stocks, but also to any other position where aspiration is required. This operation is possible and equipment is not damaged, but the aspiration power can be subjected to changes. Diameters for air passage are from 30 to 40 mm on stocks; polishing and sand-blaster machines, and retouch boxes have diameters of about 80 mm; moreover, they require a lower head. The unit (**item 090081**) complete with aspiration unit and filter (fig. 11) is appropriate for aspiration outlets with 80mm Ø. It is noiseless, 67.5 dB (A), provides an air flow of 270 cum/h and uses up 200 W.

Safety in working places

When machines are used there is always the problem of regulations and safety in working places. As already mentioned, only small or middle-sized machines can be placed in the same working room, and they must not produce annoying noises nor exhaust polluted air. In no way they must be hazardous for the health and safety of operators.

Machines must also be fitted with all protective devices, manufactured according to safety regulations, to prevent any accident caused by movements, temperature or electrical supply.

Some advice about the functional and regulation-complying features of a technological room for a dental laboratory can be find below.

Control panels and technological room in a dental laboratory

It is advisable to place aspiration units, compressors, electrical panels controlling machine operations in a duly equipped technological room.

Keep the machine room clear of anything that is not related to the machines contained therein, taking special attention of flammable materials; make sure that there is no risk of formation of corrosive, flammable or explosive mixes. Install a smoke detector for fire prevention, connected to a non-stop surveillance service.

Make sure that in this room temperature stays between a minimum of +5 °C and a max. of + 35 °C. The electric supply must be sufficient for the loads specified on machine plates and three phases + neutral + ground must be in place.

Grounding is an indispensable condition also for the single-phase. Supply lines, electric units and machines must be protected from direct and indirect electric contacts, overloads and overcurrents, in compliance with C.E.I. 64-8

regulations (available upon request) for first class equipment. It is advisable to protect the supply network from extra-currents.

Outside the room, install a fire extinguisher of appropriate size and periodically check its efficiency. Allow access to personnel only, ensure that no repair or inspection is carried out on machines in operation or ones connected to the mains. Machine casing is not only a sound-reducing device but it is also a protection from accidental contacts.

Machine installation must be carried out by an expert engineer, duly trained and equipped with the necessary tools. The engineer installing the equipment must refer to the machine manual to carry out the machine final test and to instruct users on use and routine maintenance. Technical maintenance must be carried out by engineers duly trained and equipped with the necessary tools.

Arrange for a periodical inspection of equipment and systems. This inspection is not only a means to prevent the laboratory from forced stops and to check its operating conditions, but also a way of preventing accidents and fires (according to Italian law 626).

Maintenance

Extraordinary maintenance must be carried out by engineers duly trained and specifically expert in the equipment. Maintenance operations will be scheduled considering the age and working pace of the machines; it is a good habit to keep a record of maintenance operations.

Standard maintenance can be carried out by the laboratory personnel, duly trained. Besides referring to the machine manual, which must be delivered together with the machine or with the system, the user can address to the retailer, to the installation engineer, or directly to the manufacturer.

Dust container must be emptied periodically, the personnel will schedule the operation according to experience made on first checks. Schedule must be updated according to working pace.

Filter components assembled above the paper filter must be cleaned about every three months; motorised electrovalve cover must be taken off to clean inside with an air blow twice a year.

Once a year, piping must be checked for any air leaks; after checking, start the system and make sure that aspiration is normal and that the filter automatic cleaning starts after every aspiration stop.

Important notice

- The manufacturer is at your disposal to supply spare parts, literature, instructions and anything that might be useful.
- Our Concessionaires, Agents, Retailers and authorised Service Engineers are provided with split up drawings, electrical diagrams, instructions and updating about maintenance and service.
- The guarantee and manufacturer's liability cease in case the appliances and/or the plants have been tampered with operations of any kind done by persons unfit and so unauthorised by the manufacturer.
- For any usage not contemplated or specified in this manual, please contact the manufacturer.

Transport and storage

- *Packed equipment can be transported and stored at a temperature range of -10 °C +60 °C.*
- *Packages must be kept away from water and splashing and cannot tolerate humidity >70%.*
- *Packages with the same weight can be stored in piles of three only.*

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

**CONNECTION FOR LABORATORY CENTRAL
PLANT WITH MOTORIZED VALVE**

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

DUST SEPARATOR WITH PNEUMATIC AND AUTOMATIC CLEANING OF THE FILTER

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Portata 270 m³/h.

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L= 670 P= 500 H= 830.

DUST SEPARATOR WITH PNEUMATIC AND AUTOMATIC CLEANING OF THE FILTER + ASPIRATION UNIT

Recommended for:

aspiration for polishing machine, sand blaster and box for retouch.

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L= 650 P= 700 H= 1680.

DUST SEPARATOR 75 I - MANUAL SHAKING OF THE FILTER

The filter is cleaned by a lever located outside.

Dust is collected inside a disposable plastic bag which is easy to replace.

Weight: 60 Kg.

Air flow: from 500 to 2000 m³/h.

From 6 to 24 working places in aspiration.

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L= 650 P= 700 H= 1950.

DUST SEPARATOR 75 I - PNEUMATIC CLEANING OF THE FILTER

The filter is cleaned automatically with compressed air.

Dust is collected inside a disposable plastic bag which is easy to replace.

Weight: 82 Kg.

Air flow: from 500 to 2000 m³/h.

From 6 to 24 working places in aspiration.

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ASPIRATION ASSEMBLY

(can be customized according to the client's requirements)

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
The motor is protected by a thermal device.

Output power single-phase : 0.4 kW - 3.1 A, maximum flow 1250 l/min, maximum operating head for continuous service 1300 mm H₂O. Single work place.

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

Output power single-phase: 0.75 kW - 5.5 A, output power three-phase: 0.75 kW - 3.4/1.95 A, maximum flow 2000 l/min, maximum operating head for continuous service 1400 mm H₂O. For one work station at 8-10 meters. Sound pressure level with canalized air dB (A) 67. W= 302.5 D= 302 H= 319.

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Output power single-phase: 1.1 kW - 7.6 A, output power three-phase: 1.5 kW - 6.4/3.7 A, maximum flow 3300 l/min, maximum operating head for continuous service 1500 mm H₂O. For 2 work places at 8-10 metres. Sound pressure level with canalized air 68.5 dB (A). W= 334 D= 405 H= 341.

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
Output power three-phase: 2.2 kW - 9/5.2 A, maximum flow 5000 l/min, maximum operating head for continuous service 1750 mm H₂O. For 3 work places working at the same time at 10-12 metres. Sound pressure level with canalized air 71 dB (A).

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
L= 422 P= 566 H= 443.

Output power three-phase: 4 kW - 15.6/9 A, maximum flow 8333 l/min, maximum operating head for continuous service 1700 mm H₂O.

For 4 work places working at the same time at 10-12 metres.

Sound pressure level with canalized air 76.5 dB (A).

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
L= 530 P= 530 H= 585.

Output power three-phase: 3 kW - 12/7 A, maximum flow 9500 l/min, maximum operating head for continuous service 900 mm H₂O.

For 6 work places working at the same time with the unit 20 mt away.

Sound pressure level with canalized air 78 dB (A).

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L= 570 P= 645 H= 625.


Output power three-phase: 4 kW - 15.6/9 A, maximum flow 20000 l/min, maximum operating head for continuous service 1200 mm H₂O.

For 12 work places working at the same time with the unit 25 mt away.

Sound pressure level with canalized air 78 dB (A).

W= 570 D= 645 H= 625.

Fig. 15
Draw 15



Fig. 16
Draw 16

ITALIAN PATENTS OR PATENT APPLICATIONS:

1201707 - 1234267 - 1234828 - 1259318 - 1.187.187 - 1253460 - 0766008 - 1236271 - 01242921

FOREIGN PATENTS OR PATENT APPLICATIONS:

EP 0040181 - AU 546.143 - US 4,386,910 - EP 0 638 295 - EP 0 254 687 - AU 590433 - US 4,787,846 - US 5,039,405
- EP 0335061 - US 5,002,486 - EP 0211808 - AU 580839 - US 4,684,345 - EP 0 557 251 - US 5,330,641 - EP 0766008
- US 4,710,209

PENDING PATENT

IT M093U000019 - EP 0 749 728 - IT M095U000030 - JP 168553/97 - IT M097A000139 - IT M098A000019 - IT M098A000119

CATTANI S.p.A.

Via G. Natta, 6/A - 43100 Parma (Italy)

Tel. +39 0521 607604 - Telefax +39 0521 607628-607855

<http://www.cattani.it> Email: cattani@tin.it

ESAM S.p.A.

Via G. Natta, 4/A - 43100 Parma (Italy)

Tel. +39 0521 607613 - Telefax +39 0521 607628-607855

<http://www.esam.it> Email: esamsa@tin.it